

► *Dental Stones:*

	Consistency ml/100g	Set Time	Dry Comp Strength	Setting Expansion	ISO Type
--	------------------------	-------------	----------------------	----------------------	-------------

Lab Stone - Blue, White, Yellow	30	8 - 10	8,000	0.080%	3
Lab Stone - Yellow SLOW SET	30	10 - 12	8,000	0.080%	3
High Strength - Yellow, White, Blue	28	8 - 10	9,500	0.080%	3
Orthodontic Super Stone	26	9 - 11	11,000	0.140%	3

► *Die Stones:*

Pink/White	24	10 - 12	12,500	0.067%	4
White Die Stone 16,000 psi	21	12 - 15	16,000	0.152%	5
Special Pink	21	12 - 15	16,000	0.152%	5
Green Die stone	21	12 - 15	16,000	0.180%	5
Colorburst ® Green	21	12 - 15	16,000	0.152%	5
Colorburst ® Blue	22	12 - 15	16,000	0.075%	4
Prep Stone™ Ivory	21	14 - 16	17,000	0.075%	4
Greystone XH	21	14 - 16	18,000	0.143%	5
Colorburst ® Blue XH	21	14 - 16	18,000	0.143%	5
Tan Die Stone	21	14 - 16	18,000	0.075%	4

Resinous Die Stones:

Salmon Resin	21	14 - 16	16,000	0.143%	5
Green Resin	21	14 - 16	16,000	0.143%	5
Blue Resin	21	14 - 16	16,000	0.075%	4
Ivory Resin	21	14 - 16	14,000	0.075%	4
Ivory Resin XH/ Calif. Tan	21	14 - 16	18,000	0.075%	4
Micro Resin	21	14 - 16	14,000	0.150%	5
Pacific Blue Resin	21	14 - 16	18,000	0.150%	5

► *Plasters:*

Laboratory - Regular Set	47	6 - 8	3,600	0.090%	2
Fast Set	47	3 - 4	3,600	0.090%	2
Slow Set	47	8 - 10	3,600	0.090%	2
Orthodontic Plaster	38	12 - 14	6,500	0.095%	2
Orthodontic Blend	35	14 - 16	7,200	0.110%	2
Mounting Plaster	50	1.5	5,000	0.077%	1

► *Specialty Stones:*

A & M Plaster (Flasking Stone)	38	8 - 10	6,500	0.080%	1
Quickstone (White)	21	4 - 5	14,000	0.080%	4
Optic Scanning Stone	21	14 - 16	18,000	0.075%	4
Mounting Stone	28	4 - 5	8,500	0.080%	3

Fluid Base Stones:

White, Buff & Blue	23	8 - 10	12,500	0.080%	4
--------------------	----	--------	--------	--------	---