



Premier Authorized Distributor

**ETI DIGITAL  
TECHNOLOGY**

*The Future in Dental Healthcare*

**exocad**



**UPGRADE  
NOW!**

AT A GLANCE   ENGLISH

*DentalCAD 3.3 Chemnitz*

**The most intelligent  
*DentalCAD* yet!**

DentalCAD® 3.3 Chemnitz breaks new ground in workflow efficiency and has all the tools you've been waiting for.

### AI Design\* now available for multiple units


Save time and clicks on:

- Up to three adjacent posterior crowns
- Posterior bridges with up to three elements




### Implant Module gets a major upgrade

- Better abutment proposals, auto-optimized for the respective tooth position
- Easier editing and customization
- Save and reuse custom abutment shape presets
- 3D preview of each implant component in library




### New split denture workflow, available in the Bar Module

- The new split denture workflow allows users to **design both a bar and its suprastructure** in one swift workflow
- The software splits the design into two precisely fitting components, **saving time and minimizing manual adjustments** for All-On-X hybrids



### Improved AI-based automation in bite splint design

- **Speed up** your bite splint design thanks to AI automation




## Enhanced AI for smile design

- **Drive case acceptance** with photorealistic AI smile previews with *TruSmile™ Photo\** and *TruSmile™ Video\** based on your smile design
- Automatic alignment of the 3D jaw scan and the patient smile image
- Automatic lip line detection




## Instant Anatomic Morphing now in chain mode

- **Accelerate workflows** for bridges and large tooth setups
- Real-time, synchronized adaptation of occlusal morphology for multiple restorations




## Ivoclar's clinical process integration

- Ivoclar's clinical process for digital dentures is now supported by *DentalCAD*, aiming to **reduce the number of required patient visits**
- The Ivoclar-specific workflow for functional impression-taking and bite registration utilizes the Gnathometer CAD and UTS CAD



## Integration with the DS Core cloud platform

- **Receive cases** from Dentsply Sirona and other scanners



## Typeahead search for menus

- **Find features instantly**—just start typing and jump straight to the function you need

Type to search menu options...

 Margin Line Detection

 Choose Library Teeth

 Tooth Model Placement

 Free-Forming

 Delete Reconstructions

 Adjust Scan Data Orientation

 Adapt to Pre-Op Scan

 Correct Jaw Scan Relation

 Scan Data Editor

 Correct Pre-Op Scan Placement

 3D Data Editor

## Symmetric tooth shape editing for esthetic design

- **Create symmetrical smiles efficiently** with harmonious tooth shapes





## Artex Print&Click new in the Model Creator Module

- Lets you **mount printed models** on Artex articulators without tools or plaster

## Additional highlights

- The comprehensive mesh editing tool now **supports triangle-based or clean cuts** and facilitates **gap bridging** within the mesh
- New **Inspira™ tooth shapes** for crown and bridge workflows
- New **KaVo PROTARevo 7** articulator added, bringing total number of supported articulators to 20



Premier Authorized Distributor

**ETI DIGITAL  
TECHNOLOGY**  
The Future in Dental Healthcare

**Direct 714-238-1490**

**[www.etidigitaltechnology.com](http://www.etidigitaltechnology.com)**

**[sales@etidigitaltechnology.com](mailto:sales@etidigitaltechnology.com)**



**FOLLOW US!**

**COMPLETE SUPPORT FROM OUR EXPERIENCED DIGITAL TEAM**